


# Cannulated Screw Product Series

CE 0197 ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS  
INSTRUMENT

## Headless Cannulated Compression Screw System



Product Name	Material	Diameter	Length	Product No.
Headless Cannulated Compression Screw	Titanium Alloy	Φ2.5	8、 10、 12、 14、 16、 18、 20、 22、 24、 26、 28、 30	T5702508- T5702530
		Φ3.5	20、 22、 24、 26、 28、 30、 32、 34、 36、 38、 40	T5703520- T5703540
		Φ4.5	30、 32、 34、 36、 38、 40、 42、 44、 46、 48、 50、 52、 54、 56、 58、 60	T5704530- T5704560
		Φ6.5	40、 45、 50、 55、 60、 65、 70、 75、 80、 85、 90、 95、 100、 105、 110、 115、 120	33504040- 33504120

- Hollow design, improve the accuracy of the percutaneous implant.
- Variable pitch thread and taper design pitch from big to small, the speed of the head into the bone quicker than the speed of the tail, when the screw is implanted, the fracture surface is well pressured.
- Full threaded design, better keep the pressure function when receiving periodic load.
- Headless design, to reduce collisions and stimulation to the soft tissue in the joint area in use.

## Counter Sink Cannulated Compression Screw System

**Φ2.5 Counter Sink Cannulated Compression Screw  
33507**



Product Name	Material	Diameter	Length	Product No.
Counter Sink Cannulated Compression Screw	Titanium Alloy	Φ 2.5	10	33507010
			12	33507012
			14	33507014
			16	33507016
			18	33507018
			20	33507020
			22	33507022
			24	33507024
			26	33507026
			28	33507028
30	33507030			

**Φ3.0 Counter Sink Cannulated Compression Screw  
T5503**



Product Name	Material	Diameter	Length	Product No.
Counter Sink Cannulated Compression Screw	Titanium Alloy	Φ 3.0	14	T5503014
			16	T5503016
			18	T5503018
			20	T5503020
			22	T5503022
			24	T5503024
			26	T5503026
			28	T5503028
			30	T5503030
			32	T5503032
			34	T5503034
			36	T5503036
			38	T5503038
			40	T5503040

**Φ4.0 Counter Sink Cannulated Compression Screw  
T5504**



Product Name	Material	Diameter	Length	Product No.
Counter Sink Cannulated Compression Screw	Titanium Alloy	Φ 4.0	26	T5504026
			28	T5504028
			30	T5504030
			32	T5504032
			34	T5504034
			36	T5504036
			38	T5504038
			40	T5504040
			42	T5504042
			44	T5504044
			46	T5504046
			48	T5504048
			50	T5504050

## Cannulated Compression Screw Series

**Cannulated Compression Screw(Φ3.5)  
T551**



Product Name	Material	Diameter	Length	Product No.
Cannulated Compression Screw	Titanium Alloy	Φ3.5	20、22、24、 26、28、30、 32、34、36、 38、40	T5513520- T5513540

**Cannulated Compression Screw(Φ4.0)  
T552**



Product Name	Material	Diameter	Length	Product No.
Cannulated Compression Screw	Titanium Alloy	Φ4.0	26、28、30、 32、34、36、 38、40、42、 44、46、48、 50	T5524026- T5524050

**Cannulated Compression Screw(Φ4.5)  
T553**



Product Name	Material	Diameter	Length	Product No.
Cannulated Compression Screw	Titanium Alloy	Φ4.5	24、26、28、 30、32、34、 36、38、40、 42、44、46、 48、50、52、 54、56、58、 60	T5534524- T5534560

**Femoral Neck Cannulated Compression Screw(Φ6.5)  
T562**



Product Name	Material	Diameter	Length	Product No.
Femoral Neck Cannulated Compression Screw	Titanium Alloy	Φ6.5	60、65、70、 75、80、85、 90、95、100、 105	T56265060- T56265105

Note: apply to the patient of smaller femur

**Femoral Neck Cannulated Compression Screw(Φ7.3)  
T563**



Product Name	Material	Diameter	Length	Product No.
Femoral Neck Cannulated Compression Screw	Titanium Alloy	Φ7.3	60、65、70、 75、80、85、 90、95、100、 105、110、115	T56373060- T56373115


**Dual Ends Cannulated Compression Screw(Φ6.5)  
T561**



Product Name	Material	Diameter	Length	Product No.
Dual Ends Cannulated Compression Screw	Titanium Alloy	Φ6.5	65、70、75、 80、85、90、 95、100、105、 110、115	T56165065- T56165115



## Φ2.0、Φ2.5 Self-tapping & Drilling Screw System



Product Name	Material	Diameter	Length	Product No.
Self-tapping & Drilling Screw	Titanium Alloy	Φ2.0	11、12、13、 14、15、16、 18、20	T5012011- T5012020
		Φ2.5	11、12、13、 14、15、16、 18、20、22	T5012511- T5012522

Usable range:  
Fracture of the maniphalanx  
Avulsion of trauma  
Fracture of phalange