


Kirschner Wire (Sterile Packaging)

CE 0197 ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

Intramedullary Pins (Sterile Packaging)



Character Of The System

- National class III registration certificate, more safe to implant.
- Titanium alloy material, good biocompatibility.
- Sharp triangle tip, more convenient to implant.

Metal Skeletal Pins(Double End Tips) T896



Product Name	Packaging From	Material	Diameter	Length	Product No.
Metal Skeletal Pins (Double End Tips)	Sterile Packaging	Titanium Alloy	0.8	150	T89608150-S
			1.0	150	T89610150-S
			1.2	150	T89612150-S
			1.5	150	T89615150-S
			1.8	150	T89618150-S
			2.0	150	T89620150-S
2.5	150	T89625150-S			

Metal Skeletal Pins T897



Product Name	Packaging From	Material	Diameter	Length	Product No.
Metal Skeletal Pins	Sterile Packaging	Titanium Alloy	2.0	250	T89720250-S
			2.5	250	T89725250-S
			3.0	250	T89730250-S
			3.5	250	T89735250-S
			4.0	250	T89740250-S

Sterile Packaging

